

CARTILAGE PIERCING EARRING GUIDE

Whether you are curious about our varying post lengths or how to create your ideal side profile, this guide strives to answer all of these curiosities. We recommend you to firstly read through the information in this guide, then feel free to email info@sarahandsebastian.com for any other enquiries.

Our earrings

Our internally threaded Cartilage posts are made from 14k yellow or white gold, with a 10k yellow or white gold pendant. The cartilage gauge is 18G and the backings are available in 5mm, 6,5mm and 8mm lengths.

Fitting guide

Unscrew the flat backing from the screw-in pendant by turning it anti-clockwise until the post loosens & becomes detached from the pendant. Then, insert the post into your piercing from behind & place the pendant in front of the piercing. Once the pendant is touching the post & is in position, screw the pendant clockwise into the post to ensure it is secure & tight.



The pendent with the post screwed in.




The pendent and post screwed undone.

Size guide

We believe that post lengths are a personal preference and based on each individual, however, we do have recommendations for its positioning on your ear which you can view in the image to the right.

Please note that this is a guide only and based on healed piercings on an average anatomy ear. If you have any concerns about sizing or if the piercing size is suitable for your piercing we recommend consulting with a professional beforehand.

Forward Helix	5mm
Tragus	6.5mm
Lobe	5mm-8mm
Conch	6.5mm-8mm
Flat	6.5mm-8mm



All jewellery suggestions on this page are generic. These are based on the average size of fully healed piercings on average sized anatomy. Any questions or concerns should be addressed by a professional so that you can get the ideal fit for your ears.